

# MARLEY

---

## REHAB BUGGY



marley

# REHABILITATION BUGGY



The Marley is an eye-catching rehabilitation buggy that combines simplicity and effectiveness into one quality design. The lightweight and easy-to-fold frame has been developed to meet the postural needs of each child, making the Marley a practical and functional buggy that fits easily into the daily life of any busy family.

The Marley comes with premium machine-washable black and grey upholstery, available in two sizes for excellent growth within the system and years of continued use. Made exclusively by RMS - the Marley is compatible with a range of high-quality RMS accessories and postural aids that compliment the Marley's standard features.

---

## KEY FEATURES

### STANDARD FRAME

- Available in two sizes
- Fully adjustable buggy frame
- 15° to 50° of tilt in space
- Light, compact and easy to fold
- Pelvic belt fitted as standard
- Height adjustable pram handle
- Shopping basket fitted as standard
- Full range of accessories available
- Crash tested to ISO7176-19
- Maximum user weight of 32 kg
- Quick release rear wheels
- Solid PU tyres
- Machine washable upholstery

### ACCESSORIES

- Leg abductor and hip guides
- Comfortable harnesses
- Lateral thoracic supports
- Occipital headrest
- Backrest extension
- Contour tray
- Padded grab bar
- Anklesures and foot straps
- Sun and rain covers

Build you own buggy using our new interactive product builder at [rms-rehab.co.uk](https://www.rms-rehab.co.uk)





## ERGONOMIC PRAM HANDLE

A height adjustable and comfortable handle that can be angled and moved easily with a pair of buttons. The Marley pram handle also folds down into a compact position for when the buggy isn't in use to save space.

---



## TILT IN SPACE

The Marley has a continuous tilt in space range of 15° to 50° to meet the individual requirements and comfort of the user. The tilt mechanism can be activated simply and smoothly without the need for any extra tools.

---



## SEAT BASE & PELVIC SUPPORTS

Every Marley comes with pelvic supports as standard, providing excellent postural positioning while the user is in the buggy. The seat base is also thickly padded to keep the user perfectly comfortable at all times.

---



## FOLDING FRAME

The Marley frame folds quickly and easily in one simple motion to a compact size - ideal for both transportation and storage. The rear wheels can then be removed for an even smaller storage footprint.

---



## SHOPPING BASKET

All Marley buggies come with a handy netted shopping basket as standard. This gives parents and carers a great deal of extra storage space for medical components or other items while out and about.

# TECHNICAL SPECIFICATIONS

Available in two sizes and with a wide range of adjustments - the Marley has the potential for a large range of growth - perfect for fast growing children!

	SIZE 2		SIZE 1	
	MIN	MAX	MIN	MAX
SEAT WIDTH (with pelvic supports fitted)	200 mm	270 mm	260 mm	320 mm
SEAT DEPTH	230 mm	280 mm	290 mm	350 mm
SEAT TO FLOOR HEIGHT (upright position)	-	490 mm	-	520 mm
BACKREST HEIGHT (from seat base)	-	510 mm	-	570 mm
BACKREST HEIGHT (with backrest extender)	-	625 mm	-	720 mm
THORACIC WIDTH	120 mm	210 mm	170 mm	240 mm
THORACIC HEIGHT (from seat base)	150 mm	340 mm	160 mm	430 mm
FOOTPLATE HEIGHT	130 mm	280 mm	125 mm	330 mm
FOOTPLATE DEPTH	-	160 mm	-	180 mm
PUSH HANDLE HEIGHT (from floor)	910 mm	1250 mm	925 mm	1270 mm
TILT IN SPACE ANGLE	15°	50°	15°	50°
OVERALL WIDTH	-	580 mm	-	640 mm
OVERALL DEPTH	-	790 mm	-	820 mm
OVERALL HEIGHT	-	1070 mm	-	1160 mm
BUGGY WEIGHT	15 kg	-	17 kg	-
BUGGY FOLDED SIZE	800 x 590 x 500 mm		830 x 640 x 500 mm	

## CREATE YOURS AT [WWW.RMS-REHAB.CO.UK](http://WWW.RMS-REHAB.CO.UK)

You can build your very own Marley buggy using our **NEW** interactive product builder on the RMS website. Designing your own Marley means you can get exactly what you need with all the accessories and options you require. You can check out easily and simply - getting yourself a Marley in minutes!



Thompson House, Unit 10, Styles Close,  
Sittingbourne, Kent, ME10 3BF.

**Tel:** +44 (0) 1795 477280  
**Fax:** +44 (0) 1795 229692  
**Email:** sales@rms-kent.co.uk

[www.rms-rehab.co.uk](http://www.rms-rehab.co.uk)

   @RMSRehab