



# SPRING

---

## NEWSLETTER

## NEW PRODUCT

### MY-LYTE DISCOVERY REHABILITATION BUGGY



The **NEW** My-Lyte Discovery rehab buggy incorporates the popular Shuttle Discovery seat unit onto the new My-Lyte base - a lightweight pushchair frame manufactured by RMS.

A result of extensive development - The My-Lyte frame was designed specifically for the Shuttle Discovery seat unit by our in house team of rehabilitation experts, with concepts based off our amazingly successful Marley Buggy - an industry leading postural rehab buggy.

The frame of the My-Lyte is extremely lightweight, folds in one smooth motion and will fit into the boot of most cars, thus making life easier for the parents who use it.

The modular design can be easily altered with a wide range of parts and accessories, perfect for if the child's needs evolve over time. All attachments have easy access zips so covers can be removed for simple machine washing. There is also a 'wipe-down' fabric option for messy children!

# NEW PARTNERSHIPS

## KUDU PAEDIATRIC WHEELCHAIR WITH RMS SEATING

The Kudu is manufactured by R82, so RMS have partnered up to produce a range of exclusive contoured seat cushions to offer medium to high levels of postural support and make the Kudu a more accessible option for the paediatric market.

With adjustable backrest height and seat depth - the Kudu provides excellent growth within the system, coupled with backrest and seat recline that can be adjusted easily without the need for tools.

The RMS seat cushions are available in a range of upholstery colours and a variety of contour shapes - we also supply a selection of British-made accessories, headrests, belts and harnesses that simply bolt onto the existing Kudu frame to further enhance postural support.



## CONVAID CRUISER REHABILITATION BUGGY



**The Convaid Cruiser is a rehab buggy for children requiring mild to medium levels of postural support featuring a self-tensioning seat.**

The durable, stylish and convenient Cruiser is the perfect solution for an adventurous family, folding away easily and quickly when out and about. The buggy comes complete with a 2-point padded pelvic belt; plus extendable backrest and seat depth giving this buggy up to five years of growth.

The Convaid Cruiser features a unique self-tensioning seating system which provides therapeutic pelvic positioning similar to a solid rehab chair - but all in a folding buggy. Once the user is removed, the Cruiser folds easily into a compact and lightweight package that can be easily placed into the boot of a vehicle.

**VIEW ALL OUR PRODUCTS ON THE RMS WEBSITE - [RMS-REHAB.CO.UK](http://RMS-REHAB.CO.UK)**


# NEW PARTNERSHIPS

## EASYSTAND PAEDIATRIC STANDERS

**As part of an exciting growing partnership that started at the end of 2020 - RMS have teamed up with EasyStand to supply their range of paediatric standing frames.**

These premium standers have a series of products that can cover children from 18 months old right up to large teenagers and early adulthood - along with being able to work with a range of disabilities.

### THE BANTAM



The Easystand Bantam is the ideal paediatric stander, creating a natural upright position with the head, hips, knees, and ankles all in vertical alignment. This upright position offers greater balance and is the optimal position to strengthen bones and muscles.

Another advantage of the Bantam is that it can switch quickly between posterior tilt, supine lying, hip-knee flexion and vertical standing - with a series of support pads placed across the frame to make sure the child is securely and comfortably held in every position.

A range of accessories, frame colours, and upholstery patterns; makes this stander a unique and exciting option for all children to meet their standing requirements.

### THE ZING



The Easystand Zing is the only stander on the market that allows movement from supine, through vertical and into prone positioning without having to move the tray, pads, or requiring transfer to and from the stander.

The Zing provides therapeutic opportunities not possible with other multi-position standers - use prone for head control and when fatigue sets in, simply move to a slight supine position to rest but remain in the standing position.

The Easystand Zing also comes with the capability to give up to 30 degrees of abduction to each leg, all accomplished through a single adjustment, promoting musculoskeletal development for children and infants during critical growth years.